

CAB20346

Product Information

Product SKU:	CAB20346	Gene ID:	940	Size:	20uL, 100uL
Clone No:	ARC5082-01	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	45kDa	Conjugate:	Unconjugated
Calculated MW:	25kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is essential for T-cell proliferation and survival, cytokine production, and T-helper type-2 development. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: WB,1:500 - 1:1000 FC,1:50 - 1:200

Synonyms: Tp44; CD28

Purification Method: Affinity purification

Immunogen: Recombinant protein of human CD28 Full length (P10747) .

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.