

CAB20351

Product Information

Product SKU:	CAB20351	Gene ID:	170384	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	56kDa	Conjugate:	Unconjugated
Calculated MW:	56kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable alpha-(1->3)-fucosyltransferase activity. Predicted to be involved in fucosylation. Predicted to be located in Golgi cisterna membrane. Predicted to be integral component of membrane.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: FUCTXI; FUT11

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 25-205 of human FUT11 (NP_775811.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.