

## CAB20368

---

### Product Information

<b>Product SKU:</b>	CAB20368	<b>Gene ID:</b>	4331448	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Oryza sativa

---

### Additional Information

<b>Observed MW:</b>	75kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	56kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Component of the anaphase promoting complex/cyclosome (APC/C), a cell cycle-regulated E3 ubiquitin-protein ligase complex that controls progression through mitosis and the G1 phase of the cell cycle.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** TAD1; TE; OsCCS52A; OsTE

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide of Oryza sativa OsTE.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.