

## CAB20384

---

### Product Information

|                     |          |                      |        |                    |             |
|---------------------|----------|----------------------|--------|--------------------|-------------|
| <b>Product SKU:</b> | CAB20384 | <b>Gene ID:</b>      | 242408 | <b>Size:</b>       | 20uL, 100uL |
| <b>Clone No:</b>    | -        | <b>Host Species:</b> | Rabbit | <b>Reactivity:</b> | Mouse,Rat   |

---

### Additional Information

|                       |       |                   |              |
|-----------------------|-------|-------------------|--------------|
| <b>Observed MW:</b>   | 55kDa | <b>Conjugate:</b> | Unconjugated |
| <b>Calculated MW:</b> | 54kDa | <b>Isotype:</b>   | IgG          |

---

### Immunogen Information

**Background:** Orthologous to human FAM221B (family with sequence similarity 221 member B).

**Recommended Dilution:** WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

**Synonyms:** 4930412F15Rik; Fam221b

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 280-487 of mouse Fam221b (NP\_780726.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.