

CAB20413

Product Information

Product SKU:	CAB20413	Gene ID:	89790	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	100kDa	Conjugate:	Unconjugated
Calculated MW:	77kDa	Isotype:	IgG

Immunogen Information

Background: SIGLECs are members of the immunoglobulin superfamily that are expressed on the cell surface. Most SIGLECs have 1 or more cytoplasmic immune receptor tyrosine-based inhibitory motifs, or ITIMs. SIGLECs are typically expressed on cells of the innate immune system, with the exception of the B-cell expressed SIGLEC6 (MIM 604405).

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: SLG2; PRO940; SIGLEC-10; SIGLEC10

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 627-697 of human SIGLEC10 (NP_149121.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.