

CAB20439

Product Information

Product SKU:	CAB20439	Gene ID:	401505	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	6kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be involved in protein targeting to mitochondrion. Located in mitochondrion. Part of mitochondrial outer membrane translocase complex.

Recommended Dilution: IF/ICC, 1:50 - 1:200

Synonyms: Tom5; C9orf105; bA613M10.3; TOMM5

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-51 of human TOMM5 (NP_001001790.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.