

CAB20447

Product Information

Product SKU: CAB20447	Gene ID: 9975	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Mouse,Rat

Additional Information

Observed MW: 70kDa	Conjugate: Unconjugated
Calculated MW: 65kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a member of the nuclear hormone receptor family, specifically the NR1 subfamily of receptors. The encoded protein functions as a transcriptional repressor and may play a role in circadian rhythms and carbohydrate and lipid metabolism. Alternatively spliced transcript variants have been described.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: RVR; BD73; EAR-1R; REVERBB; REVERBbeta; Rev-erb beta/NR1D2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 220-300 of human Rev-erb beta/NR1D2 (NP_005117.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.