

## CAB20455

---

### Product Information

<b>Product SKU:</b>	CAB20455	<b>Gene ID:</b>	2357	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	Refer to figures	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	38kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene encodes a G protein-coupled receptor of mammalian phagocytic cells that is a member of the G-protein coupled receptor 1 family. The protein mediates the response of phagocytic cells to invasion of the host by microorganisms and is important in host defense and inflammation.

**Recommended Dilution:** WB,1:500 - 1:2000 IHC-P,1:50 - 1:200

**Synonyms:** FPR; FMLP; FPR1

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 250-350 of human FPR1 (NP\_002020.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.