

CAB20463

Product Information

Product SKU:	CAB20463	Gene ID:	8572	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	35kDa	Conjugate:	Unconjugated
Calculated MW:	35kDa	Isotype:	IgG

Immunogen Information

Background: Enables alpha-actinin binding activity; protein homodimerization activity; and protein phosphatase binding activity. Involved in actin cytoskeleton reorganization. Located in several cellular components, including lamellipodium; perinuclear region of cytoplasm; and stress fiber. Part of filamentous actin. Implicated in osteoporosis.

Recommended Dilution: WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200

Synonyms: RIL; PDLIM4

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 80-210 of human PDLIM4 (NP_003678.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.