

CAB20472

Product Information

Product SKU:	CAB20472	Gene ID:	15894	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	85kDa	Conjugate:	Unconjugated
Calculated MW:	59kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes an integral membrane protein that binds leukocyte adhesion protein LFA-1. It participates in the innate immune response.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: CD54; Ly-47; Icam-1; MALA-2; ICAM1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 120-320 of mouse ICAM1 (NP_034623.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.