

CAB20476

Product Information

Product SKU:	CAB20476	Gene ID:	23295	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	70kDa	Conjugate:	Unconjugated
Calculated MW:	61kDa	Isotype:	IgG

Immunogen Information

Background: Enables ubiquitin-protein transferase activity. Involved in endosome to lysosome transport; negative regulation of signal transduction; and protein monoubiquitination. Located in several cellular components, including early endosome; endoplasmic reticulum; and nucleus.

Recommended Dilution: WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200

Synonyms: RNF156; MGRN1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-538 of human MGRN1 (O60291).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.