

CAB20479

Product Information

| | | |
|------------------------------|-----------------------------|------------------------------------|
| Product SKU: CAB20479 | Gene ID: 54492 | Size: 20uL, 100uL |
| Clone No: - | Host Species: Rabbit | Reactivity: Human,Mouse,Rat |

Additional Information

| | |
|-----------------------------|--------------------------------|
| Observed MW: 59kDa | Conjugate: Unconjugated |
| Calculated MW: 59kDa | Isotype: IgG |

Immunogen Information

Background: Predicted to enable ubiquitin protein ligase activity. Predicted to be involved in ubiquitin-dependent endocytosis. Located in actin cytoskeleton and cytosol.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: neur2; NEURL3; RNF67B; hNeur2; NEURL1B

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 270-440 of human NEURL1B (NP_001136123.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.