

CAB20485

Product Information

Product SKU:	CAB20485	Gene ID:	79828	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	33kDa	Conjugate:	Unconjugated
Calculated MW:	33kDa	Isotype:	IgG

Immunogen Information

Background: Enables mRNA methyltransferase activity. Involved in mRNA methylation. Predicted to be located in cytoplasm and nucleus.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: TIP; METTL8

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human METTL8 (NP_001308089.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.