

CAB20490

Product Information

Product SKU: CAB20490	Gene ID: 91607	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human

Additional Information

Observed MW: 100kDa	Conjugate: Unconjugated
Calculated MW: 103kDa	Isotype: IgG

Immunogen Information

Background: Enables tRNA binding activity. Involved in several processes, including defense response to virus; negative regulation of G1/S transition of mitotic cell cycle; and replication fork arrest. Located in cytosol; nucleoplasm; and site of DNA damage.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: SLFN8/9; SLFN11

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 132-328 of human SLFN11 (NP_689483.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.