

## CAB20499

---

### Product Information

<b>Product SKU:</b>	CAB20499	<b>Gene ID:</b>	1840	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

---

### Additional Information

<b>Observed MW:</b>	67-72kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	67kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Studies in Drosophila have identified this gene as encoding a positive regulator of the Notch-signaling pathway. The human gene encodes a protein of unknown function; however, it may play a role in basic helix-loop-helix transcription factor activity.

**Recommended Dilution:** WB, 1:100 - 1:500

**Synonyms:** hDx-1; RNF140; DTX1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human DTX1 (NP\_004407.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.