

CAB20523

Product Information

Product SKU:	CAB20523	Gene ID:	4642	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	116kDa	Conjugate:	Unconjugated
Calculated MW:	116kDa	Isotype:	IgG

Immunogen Information

Background: Enables protein domain specific binding activity. Predicted to be involved in actin filament organization; early endosome to recycling endosome transport; and vesicle transport along actin filament. Located in extracellular exosome.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: myr4; PPP1R108; MYO1D

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 180-380 of human MYO1D (NP_001290209.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.