

## CAB20533

---

### Product Information

<b>Product SKU:</b>	CAB20533	<b>Gene ID:</b>	124599	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	23kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	23kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** CD300LB is a nonclassical activating receptor of the immunoglobulin (Ig) superfamily expressed on myeloid cells (Martinez-Barriocanal and Sayos, 2006 [PubMed 16920917]).

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** CLM7; CLM-7; IREM3; TREM5; CD300b; IREM-3; TREM-5; CMRF35-A2; CD300LB

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-201 of human CD300LB (NP\_777552.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.