

CAB20534

Product Information

Product SKU:	CAB20534	Gene ID:	283742	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	48kDa	Conjugate:	Unconjugated
Calculated MW:	46kDa	Isotype:	IgG

Immunogen Information

Background: Enables identical protein binding activity and protein methyltransferase activity. Involved in positive regulation of cell population proliferation; positive regulation of gene expression; and protein methylation. Located in cytoplasm and nucleoplasm. Part of tRNA-splicing ligase complex.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: -

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-225 of human FAM98B (NP_775882.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.