

CAB20564

Product Information

Product SKU:	CAB20564	Gene ID:	404093	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	42kDa	Conjugate:	Unconjugated
Calculated MW:	42kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable ubiquitin binding activity.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: -

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 187-386 of human CUEDC1 (NP_001258804.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.