

## CAB20572

---

### Product Information

<b>Product SKU:</b>	CAB20572	<b>Gene ID:</b>	4341678	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Oryza sativa

---

### Additional Information

<b>Observed MW:</b>	52kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	44kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** This gene is mainly expressed in roots before flowering (RPKM288.3) and leaves after flowering (RPKM21.0).Modulates plant innate immunity against X.oryzae pv. oryzae (Xoo) (PMID: 21961049).Regulates negatively the basal defense responses to the compatible fungus M.oryzae (PMID: 23462973, PMID: 21726398).

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** OsWRKY28

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of Oryza sativa OsWRKY28.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.