

CAB20573

Product Information

Product SKU:	CAB20573	Gene ID:	4340170	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	67kDa	Conjugate:	Unconjugated
Calculated MW:	47kDa	Isotype:	IgG

Immunogen Information

Background: Trans-acting factor that binds specifically to the consensus nucleotide sequence 5'-TNCGTACAA-3' (By similarity). May be involved in panicle development.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: MPK1; sipk; MAPK6; OsSIPK; OsMAPK6; OsSPL16

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa OsSPL16.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.