

CAB20574

Product Information

Product SKU: CAB20574	Gene ID: 4339440	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Oryza sativa

Additional Information

Observed MW: 70kDa	Conjugate: Unconjugated
Calculated MW: 47kDa	Isotype: IgG

Immunogen Information

Background: contributes to the modification of pre-tRNA^{His} and is required for normal growth under high-temperature conditions in rice.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: AET1; tRNA^{His} guanylyltransferase

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa AET1.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.