

CAB20575

Product Information

Product SKU:	CAB20575	Gene ID:	4334553	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	35kDa	Conjugate:	Unconjugated
Calculated MW:	35kDa	Isotype:	IgG

Immunogen Information

Background: Transcription factor that possesses transactivation activity . Transcription activator involved in response to abiotic stresses. Plays a positive role during dehydration and salt stress. Binds specifically to the 5'-CATGTG-3' motif found in promoters of stress-responsive genes

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: NAC; NAC19; NAC20; snac1; SNAC1

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa SNAC1.

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.