

CAB20580

Product Information

Product SKU:	CAB20580	Gene ID:	4340427	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	107kDa	Conjugate:	Unconjugated
Calculated MW:	107kDa	Isotype:	IgG

Immunogen Information

Background: Receptor kinase involved in the regulation of thermotolerance. Functions as positive regulator of heat tolerance. May be involved in the regulation of cell proliferation and cell growth

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: OsERECTA; OsER1

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa OsERECTA.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.