

CAB20583

Product Information

Product SKU:	CAB20583	Gene ID:	4345634	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	29kDa	Conjugate:	Unconjugated
Calculated MW:	71kDa	Isotype:	IgG

Immunogen Information

Background: Is associated with a DNA binding complex that binds to the G box, a well-characterized cis-acting DNA regulatory element found in plant genes.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: GF14C; OJ1124_B05.7; GF14-C

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa GF14-C

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.