

CAB20586

Product Information

Product SKU:	CAB20586	Gene ID:	4337998	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	25kDa	Conjugate:	Unconjugated
Calculated MW:	21kDa	Isotype:	IgG

Immunogen Information

Background: WRKY genes encode transcription factors that are involved in the regulation of various biological processes

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: OsWRKY67

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa OsWRKY67

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.