

CAB20587

Product Information

Product SKU: CAB20587	Gene ID: 4341978	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Oryza sativa

Additional Information

Observed MW: 110kDa	Conjugate: Unconjugated
Calculated MW: 128kDa	Isotype: IgG

Immunogen Information

Background: Auxin response factors (ARFs) are transcriptional factors that bind specifically to the DNA sequence 5'-TGCTC-3' found in the auxin-responsive promoter elements (AuxREs).

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: ARF19; ARF7A; OsARF7a; OJ1215_E11.1; OsARF19

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa OsARF19

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.