

## CAB20594

---

### Product Information

<b>Product SKU:</b>	CAB20594	<b>Gene ID:</b>	14254601	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	MERS-CoV

---

### Additional Information

<b>Observed MW:</b>	47kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	45kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Packages the positive strand viral genome RNA into a helical ribonucleocapsid (RNP) and plays a fundamental role during virion assembly through its interactions with the viral genome and membrane protein M. Plays an important role in enhancing the efficiency of subgenomic viral RNA transcription as well as viral replication.

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** -

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-413 of coronavirus Nucleoprotein (YP\_009047211.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.