

CAB20595

Product Information

Product SKU:	CAB20595	Gene ID:	57217	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	96kDa	Conjugate:	Unconjugated
Calculated MW:	96kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a protein containing tetratricopeptide repeats. Mutations in this gene disrupt intestinal development and can cause early onset inflammatory bowel disease and intestinal atresia. Alternative splicing results in multiple transcript variants.

Recommended Dilution: WB,1:100 - 1:500

Synonyms: TTC7; GIDID; MINAT; TTC7A

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TTC7A (NP_001275880.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.