

CAB20612

Product Information

Product SKU:	CAB20612	Gene ID:	3200423	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	HCoV-HKU1

Additional Information

Observed MW:	50-60kDa	Conjugate:	Unconjugated
Calculated MW:	49kDa	Isotype:	IgG

Immunogen Information

Background: Packages the positive strand viral genome RNA into a helical ribonucleocapsid (RNP) and plays a fundamental role during virion assembly through its interactions with the viral genome and membrane protein M. Plays an important role in enhancing the efficiency of subgenomic viral RNA transcription as well as viral replication.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: -

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 100-200 of coronavirus Nucleoprotein (YP_173242.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.