

CAB20613

Product Information

Product SKU:	CAB20613	Gene ID:	3200426	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	HCoV-HKU1

Additional Information

Observed MW:	75kDa/110kDa/145kDa	Conjugate:	Unconjugated
Calculated MW:	152kDa	Isotype:	IgG

Immunogen Information

Background: [Spike protein S1]: attaches the virion to the cell membrane by interacting with host receptor, initiating the infection.

Recommended Dilution: WB, 1:100 - 1:500

Synonyms: -

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 100-200 of coronavirus Spike S1 (YP_173238.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.