

CAB20616

Product Information

Product SKU:	CAB20616	Gene ID:	39105218	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	HCoV-OC43

Additional Information

Observed MW:	75kDa/145kDa	Conjugate:	Unconjugated
Calculated MW:	150kDa	Isotype:	IgG

Immunogen Information

Background: S1 attaches the virion to the cell membrane by interacting with sialic acid-containing cell receptors, initiating the infection.; [Spike protein S1]: attaches the virion to the cell membrane by interacting with host receptor, initiating the infection.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: -

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 782-881 of coronavirus Spike S2 (YP_009555241.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.