

CAB20622

Product Information

Product SKU:	CAB20622	Gene ID:	4324637	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	43kDa	Conjugate:	Unconjugated
Calculated MW:	43kDa	Isotype:	IgG

Immunogen Information

Background: encodes an enoyl-CoA hydratase/isomerase, increases the grain yield of rice by enhancing grain number per panicle without a negative effect on the number of panicles per plant or grain weight.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: NOG1

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa NOG1.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.