

CAB20624

Product Information

Product SKU:	CAB20624	Gene ID:	4341780	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	24kDa	Conjugate:	Unconjugated
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

Background: encodes a plant-specific plant AT-rich sequence- and zinc-binding transcription factor that regulates rice grain length and spikelet number.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: GL6; SG6

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa GL6.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.