

CAB20625

Product Information

Product SKU:	CAB20625	Gene ID:	4336066	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	30kDa	Conjugate:	Unconjugated
Calculated MW:	30kDa	Isotype:	IgG

Immunogen Information

Background: Is associated with a DNA binding complex that binds to the G box, a well-characterized cis-acting DNA regulatory element found in plant genes.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: GF14B; B0812A04.1; GF14b

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa GF14b.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.