

CAB20627

Product Information

| | | | | | |
|---------------------|----------|----------------------|---------|--------------------|--------------|
| Product SKU: | CAB20627 | Gene ID: | 4340462 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Oryza sativa |

Additional Information

| | | | |
|-----------------------|-------|-------------------|--------------|
| Observed MW: | 50kDa | Conjugate: | Unconjugated |
| Calculated MW: | 35kDa | Isotype: | IgG |

Immunogen Information

Background: Transcription factor involved in abscisic acid (ABA) signaling pathway Transcription factor activity is fully activated by ABA.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: P0701E03.2; bZIP46

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa bZIP46.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.