

CAB20629

Product Information

Product SKU:	CAB20629	Gene ID:	4330221	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	51kDa	Conjugate:	Unconjugated
Calculated MW:	41kDa	Isotype:	IgG

Immunogen Information

Background: Removes the phosphate from trehalose 6-phosphate to produce free trehalose. Trehalose accumulation in plant improves abiotic stress tolerance.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: TPP1; OsTPP1; OsJ_07820; P0516F12.9

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa OsTPP1.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.