

CAB20633

Product Information

Product SKU:	CAB20633	Gene ID:	4334396	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	70kDa	Conjugate:	Unconjugated
Calculated MW:	79kDa	Isotype:	IgG

Immunogen Information

Background: EL1 is a novel and key regulator of the GA response, can specifically phosphorylates the rice DELLA protein SLR1

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: CKI; EL1; Hd16; OsJ_12923

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa Hd16.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.