

CAB20634

Product Information

Product SKU: CAB20634	Gene ID: 4335151	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Oryza sativa

Additional Information

Observed MW: 90kDa	Conjugate: Unconjugated
Calculated MW: 90kDa	Isotype: IgG

Immunogen Information

Background: Bph3 is a cluster of three genes encoding plasma membrane-localized lectin receptor kinases (OsLecRK1-OsLecRK3)

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: OsLecRK3; Bph3

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa Bph3.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.