

CAB20636

Product Information

| | | | | | |
|---------------------|----------|----------------------|---------|--------------------|--------------|
| Product SKU: | CAB20636 | Gene ID: | 4350610 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Oryza sativa |

Additional Information

| | | | |
|-----------------------|--------|-------------------|--------------|
| Observed MW: | 110kDa | Conjugate: | Unconjugated |
| Calculated MW: | 94kDa | Isotype: | IgG |

Immunogen Information

Background: Auxin response factors (ARFs) are transcriptional factors that bind specifically to the DNA sequence 5'-TGCTC-3' found in the auxin-responsive promoter elements (AuxREs).

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: arf1; ARF23; OsARF1; OsARF23

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa OsARF23.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.