

## CAB20640

---

### Product Information

<b>Product SKU:</b> CAB20640	<b>Gene ID:</b> 4335004	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Oryza sativa

---

### Additional Information

<b>Observed MW:</b> 120kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 120kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** CTB4a (cold tolerance at booting stage), encoding a conserved leucine-rich repeat receptor-like kinase. The different CTB4a alleles confer distinct levels of cold tolerance and selection for variation in the CTB4a promoter region has occurred on the basis of environmental temperature.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** LRR-PK1; CTB4a

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of Oryza sativa CTB4a.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.