

CAB20642

Product Information

Product SKU:	CAB20642	Gene ID:	4345993	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	103kDa	Conjugate:	Unconjugated
Calculated MW:	103kDa	Isotype:	IgG

Immunogen Information

Background: Plant lipoxygenase may be involved in a number of diverse aspects of plant physiology including growth and development, pest resistance, and senescence or responses to wounding.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: LOX8; LOX2.1; CM-LOX1; lox2osPil; OsHI-LOX

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa OsHI-LOX.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.