

CAB20646

Product Information

Product SKU:	CAB20646	Gene ID:	25925	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Rat

Additional Information

Observed MW:	148kDa	Conjugate:	Unconjugated
Calculated MW:	148kDa	Isotype:	IgG

Immunogen Information

Background: Enables protein domain specific binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to act upstream of or within neuron fate commitment. Located in nucleus.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: EHZF; Evi3; ZNF521

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 210-420 of human ZNF521 (NP_056276.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.