

CAB20650

Product Information

Product SKU:	CAB20650	Gene ID:	113791	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	43kDa	Conjugate:	Unconjugated
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

Background: Enables phosphatidylinositol 3-kinase catalytic subunit binding activity. Involved in negative regulation of phosphatidylinositol 3-kinase activity and negative regulation of phosphatidylinositol 3-kinase signaling. Predicted to be located in plasma membrane. Predicted to be integral component of membrane.

Recommended Dilution: WB, 1:100 - 1:500

Synonyms: HGFL; TrIP; hHGFL(S); PIK3IP1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 24-160 of human PIK3IP1 (NP_443112.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.