

## CAB20655

---

### Product Information

<b>Product SKU:</b> CAB20655	<b>Gene ID:</b> -	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Oryza sativa

---

### Additional Information

<b>Observed MW:</b> 28kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 28kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** An-1 positively regulates awn elongation, but negatively regulates grain number per panicle. Genetic variations in the An-1 locus were found to be associated with awn loss in cultivated rice.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** An-1

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of Oryza sativa An-1.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.