

## CAB20658

---

### Product Information

<b>Product SKU:</b>	CAB20658	<b>Gene ID:</b>	4330621	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Oryza sativa

---

### Additional Information

<b>Observed MW:</b>	30kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	27kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Transcriptional factor that targets the CARG motif 5'-C(A/T)TTAAAAAG-3' in the promoter of D14. Directly suppresses D14 expression to control the outgrowth of axillary buds.

**Recommended Dilution:** WB, 1:500 - 1:2000

**Synonyms:** MADS57; OsMADS57; OsJ\_08264

**Purification Method:** Affinity purification

**Immunogen:** Recombinant protein of Oryza sativa OsFD2.

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.