

CAB20661

Product Information

Product SKU: CAB20661	Gene ID: 4331983	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Oryza sativa

Additional Information

Observed MW: 36kDa	Conjugate: Unconjugated
Calculated MW: 34kDa	Isotype: IgG

Immunogen Information

Background: Involved in strigolactone (SL) signaling pathway. May function downstream of SL synthesis, as a component of hormone signaling or as an enzyme that participates in the conversion of SL to the bioactive form.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: D14; D88; HTD2

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa D14.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.