

CAB20663

Product Information

Product SKU:	CAB20663	Gene ID:	4340725	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	59kDa	Conjugate:	Unconjugated
Calculated MW:	59kDa	Isotype:	IgG

Immunogen Information

Background: a master regulator of the growth–defense trade-off in rice

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: Aldh2b; OsALDH2B1

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa OsALDH2B1.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.