

CAB20668

Product Information

Product SKU: CAB20668	Gene ID: -	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Oryza sativa

Additional Information

Observed MW: 100kDa	Conjugate: Unconjugated
Calculated MW: 48kDa	Isotype: IgG

Immunogen Information

Background: conferring resistance to BPH at seedling and maturity stages of the rice plant

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: -

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa BPH14.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.