

CAB20669

Product Information

Product SKU: CAB20669	Gene ID: 4349543	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Oryza sativa

Additional Information

Observed MW: 122.5kDa	Conjugate: Unconjugated
Calculated MW: 123kDa	Isotype: IgG

Immunogen Information

Background: Repressor of strigolactones (SL) signaling. Subjected to a negative feedback control of SL signaling.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: D53; DWARF53; OsJ_32651

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa D53.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.